

HARDENER TESTER

It comprise precisionly engineered Dial gauge, Aluminum Body and finely machined components. Dial comprises main pointer and follow pointer. Tablet is to be positioned in between with two anvils. Lower anvil stationery while upper one is moving with the help of screw exerting and unit force linear. Force is recorded on Dia Gauge by follow pointer while main pointer goes back to zero. It can be validated by dead weight upto 5 kgs, which will ensure you repeatability. It can either be fixed on bench or tableting machine.

Note : Images Shown here are illustrative. As the design & manufacturing of Machines are subject to improvement, the product supplied will be as per our Techno-Commercial offer.